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




General product information

Acryl/B are cold polymerisates based on polymethylmethacrylate (PMMA) in different colors for manufacturing BTE earmolds and ITE shells. Polymerization must be conducted under pressure in the water bath. Manual pressure pots (Aquapress REF 49800 / Acrylmat REF 53600) as well as automatic polymerization units (Polymat REF 54200) can be used to serve as pressure polymerization units. Please carefully observe the user manual for the device you are using.

Package item: See catalog

Stability/batch: See tub

Storage:

Polymer (Powder): Store the material in a dark und dry place. Storage is possible 10 years in closed container.

Monomer (Liquid): Keep cool (not over 25 °C). Storeage is possible 24 months from date of production with originally closed bottle. Ensure good ventilation in the entire room even the floor (fumes are heavier than air). Suction system at the workplace is recommended. Keep away from sources of ignition. No smoking. Protect from contamination and sunrays as otherwise premature polymerization may be caused.

Note: Close container immediately after use. Avoid direct contact with liquid material! Before processing and shaping the material, hands should be protected by wearing suitable gloves.

Processing

1. Combine monomer (liquid) and polymer (powder) to a mixing ratio of 1:2. The monomer should be cold (keep in refrigerator) for the ideal swelling behavior. We recommend mixing the material under a suction unit.

Tip: First fill the monomer into the mixing container (REF 47700 ff) and then add the polymer to avoid the formation of clusters. Use a vibrator to more easily bring up the bubbles that form while the material is mixed.

2a. Producing BTE earmolds

Different materials (plaster/duplicating silicone/duplicating material) can be used for the production of negative molds for BTE earmolds. After mixing the acrylic material to the correct ratio on the activated vibrator, gradually fill the material into the negative mold until it is completely filled. Let stand on the vibrator briefly to let possible bubbles come to the surface.

2b. Producing ITE shells, beige

Duplicating silicone transparent (REF 41000 ff) or duplicating silicone pink (REF 41801 ff) can be used for the production of negative molds for ITE shells. Warm up the duplicating silicone negative mold in the hot air oven (REF 68301) at a temperature of 60 °C (140 °F) for 15 minutes. After mixing the acrylic material to the correct ratio on the activated vibrator, gradually fill the material into the negative mold until it is completely filled. After a few seconds, a shell wall will form. Turn the filled negative mold over so the remaining liquid acrylic material can be emptied.

3. After a few seconds a matted surface will form on the acrylic material. Place the negative form into a pressure pot filled with water (40 °C/104 °F) for polymerization at 2–5 bar for 12–15 minutes.

3. After a few seconds a matted surface will form on the acrylic material. Place the negative mold into a pressure pot filled with water (40 °C/104 °F) for polymerization at 2–5 bar for 12–15 minutes.

Attention: The negative mold must be covered completely with water; otherwise, bubbles will form in the blank.

4. After polymerization remove the BTE blank or ITE shell from the negative form and shape with suitable special cutters.

5. Afterwards you can coat the earmold/ITE shell with egger LP/H lacquer (REF 30600 ff), egger LP/H lacquer antibac (REF 30621 ff) or egger LP/W lacquer (REF 30700 ff).