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GB Disinfectant wipes for hearing systems and earmolds. Effective against bacteria, viruses (HBV/HIV) and fungi (Candida albicans). Reaction time: 1 minute. Rub

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hearing systems and earmolds with amoiest disinfecting wipe (image 1) and then rub dry with an absorbent cloth after a 1 minute reaction time (image 2).

Caution: Harmful to aquatic life with long lasting effects. Contains benzisothiazolionone. May produce an allergic reaction. If medical advice is needed, have product container or label at hand. Keep out of reach of children. Read label before use. Avoid release to the environment. Do not use after the expiration date has passed.

△ **Following disinfection, always rub objects dry (e.g. hearing aids, earmolds) that come into contact with skin using an absorbent cloth.**

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

Quaternary ammonium compounds, benzy-C12-16-alkyldimethyl chloride 0.25g, Didecyldimethylammoniumchloride 0.25g, Quaternary ammonium compounds, C12-14-alkyl[(ethyphenyl)methyl]dimethyl chloride 0.25g/100g

Desinfectant wipes for hearing systems and earmolds. Effective against bacteria, viruses (HBV/HIV) and fungi (Candida albicans). Reaction time: 1 minute. Rub

Desinfectant wipes for hearing systems and earmolds. Effective against bacteria, viruses (HBV/HIV) and fungi (Candida albicans). Reaction time: 1 minute. Rub

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GB Disinfectant spray for hearing systems and earmolds. Effective against bacteria, viruses (HBV/HIV) and fungi (Candida albicans). Apply until completely moistened

(figure 1) and rub dry with an absorbent cloth after an exposure time of 1 minute (figure 2). Do not spray into the openings of the microphone and receiver! **Caution:** Harmful to aquatic life with long lasting effects. Contains benzisothiazolinone. May produce an allergic reaction. If medical advice is needed, have product container or label at hand. Keep out of reach of children. Read label before use. Avoid release to the environment. Do not use after the expiration date has passed.

△ Following disinfection, always rub objects dry (e.g. hearing aids, earmolds) that come into contact with skin using an absorbent cloth.

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